6 Steps to Successful Drenching using the **Drench-Mate**







Prepare the Drench

- Read the drench directions.
- Mix the drench powder thoroughly with warm water (fig 1).
- Close the *Drench-Mate* lid to prevent spillage (fig 2).



Prepare the Cow

- Ensure the cow is restrained (fig 1).
- Keep the cow's head low to ensure a direct path to the rumen (fig 1).
- Disengage the *Drench-Mate* nose tong and insert frick tube carefully (fig 2 and 3).
- Engage nose tong to hold tube in place (fig 4).









Figure 4

Insert Hose Correctly

- Carefully insert *Drench-Mate* hose until you feel resistance in the back of the cow's mouth. Pause and let the cow swallow (fig 1). DO NOT FORCE. Once the hose has been swallowed, feel the throat and make sure there are two hard structures. One will be the hose in the esophagus, the other will be the cow's trachea.
- Push the hose in all the way to the red mark. This will help ensure you are in the rumen (fig 2).



Distributor

Drench-Mate 9314 Swanson Rd., Sumas, WA 98295

phone: 360-988-5020 • fax: 360-988-5034









Figure 1



Listen & Smell

- Before attaching the end of the hose to the *Drench-Mate*_® pump, hold to your ear and listen (fig 1). You should hear a bubbling or gurgling sound. You should not hear breathing.
- Smell the end of the hose. You should detect the fermenting odor of the rumen (fig 2).



Figure 2



Pump Correctly



Figure 1



Figure 2

- Connect the hose to the *Drench-Mate* pump and begin pumping (fig 1). An easy, steady stroke should be used.
- With the Single Dose, tilt the bucket to the inlet and swish a few times when it is almost empty (fig 2).



Drain & Remove Hose

- Disconnect hose from pump and hold high for 10 seconds so that any remaining fluid drains out (fig 1).
- Blow any remaining fluid from the hose (fig 2).
- Disconnect from nose tong and remove hose carefully. Then remove the frick tube (fig 3).



Figure 1



Figure 2



Figure 3



Distributor Drench-Mate 9314 Swanson Rd. Sumas, WA 98295

phone: 360-988-5020 fax: 360-988-5034